

Pichia Strains

ATUM offers a selection of *Pichia pastoris* expression strains for research use. These strains are ready to use with all ATUM *Pichia* vectors. Three IP-Free *Pichia* Strains are available – wild type, slow methanol metabolizing and protease deficient strains allowing you to optimize expression in *Pichia*.

Pichia Strains

Cat #	Name
PPS-9010	wild type
PPS-9011	aox1Δ (Mut ^S)
PPS-9016	pep4Δ, prb1Δ
PPS-KT (3 strains)	Strain kit (3 strains)

200µl Package size. Strains are shipped at room temperature in YPD + 1M Sorbitol. Store at 4°C on arrival.

Protocols for strain propagation:

Strains are in YPD + 1M Sorbitol broth and are shipped at ambient temperature. When received, strains should be stored at 4°C and revived within 1 week by streaking for single colonies to YPD Agar plates (Teknova, Cat. No. Y1000 with 1% yeast extract, 2% peptone, 2% glucose and 2% agar) and incubating for 2 days at 30°C.

Master stocks are made by culturing a single colony of each strain overnight in YPD broth (Teknova, Cat No. Y5000 with 1% yeast extract, 2% peptone, 2% glucose). Harvest cells (5 min at 500 x g) and suspend cell pellet in YPD Broth containing 30% glycerol. Cells are frozen on a dry ice/ethanol bath and stored at -80C.

For Competent cell preparation, Culture and Induction protocols, see

www.atum.bio/wp-content/uploads/2013/04/Pichia_culture_induction_protocol.pdf

Intellectual Property Statement

Proteins and peptides made with *Pichia* strains (the “*Pichia* Strains”) can be used commercially without license obligations. However, neither the *Pichia* Strains, nor any strains derived from the *Pichia* Strains, may be transferred or sold to third parties, resold, modified for resale, or used to provide a service of any kind to third parties, including, without limitation, reporting the results of customer activities for a fee or other form of consideration.

Available online: <https://www.atum.bio/company/terms-and-conditions>

These strains have been jointly developed by BioGrammatics and ATUM (formerly DNA2.0).