

Bolt PCR Cloning Kit

Bolt Cloning

ATUM_{SM} (formerly DNA2.0) has developed a simple, fast, and efficient method for cloning PCR amplified products into a plasmid vector in a 5-minute bench-top reaction. This cloning kit contains an optimized cloning mix with a linearized vector that allows cloning of any PCR product amplified with proofreading DNA polymerases (not compatible with *Taq* polymerase).

Kit Components (20 Rx)

Linearized pJ201 (kan^R) in BKT-01 **or pJ204 (amp^R) vector** in BKT-02 is supplied as a 10X mix, to be diluted to 1X in the final reaction mix (10X mix contains 20ng/μl linearized vector DNA and topoisomerase).

Positive Control is a PCR amplicon of KringleYFP linked to a tet promoter, 25 ng/μl, 864 base pairs, which can be cloned in either orientation.

Sequencing Primers: 1 μM stock, use 1 μl / sequencing reaction


Forward primer: 5'- TGGTAGTGTGGGGACTC-3'

Reverse primer: 5'- TTGTCAGAATATTTAAGGGCG-3'

Storage

Store all kit components at -20°C.

Vector map



Cloning Information

To PCR your ORF: Bolt vectors do not contain multiple cloning sites (MCS), so flanking restriction sites must be incorporated in your PCR primers for easy subcloning into a vector of choice.

Note: Include 5 additional bases outside of your flanking restriction sites by adding 5 bp to each 5' end of the PCR primers. This will accommodate potential end deletions.

Bolt reaction is set up as follows:

COMPONENT	VOLUME (μl)
PCR product (~100ng)/Control	4
Linearized Bolt Vector (20ng) mix	1
Sterile ddH ₂ O	5
Total Volume	10

1. Combine components as listed above in single 1.5 ml tube. Incubate at room temperature for 5-20 minutes. For best results, pre-equilibrate PCR reaction to room temperature for 5 minutes prior to adding to vector.
2. Transform 2.5 μ l of each reaction into 50 μ l competent cells. Heat shock if appropriate, and then add 500 μ l SOC Broth. Shake for 1 hour at 37C.
3. Plate 100 μ l on LB + selection antibiotic – pJ201 on LB agar + 30 μ g/ml kanamycin; pJ204 on LB agar + 100 μ g/ml ampicillin or carbenicillin.
4. Incubate plates overnight at 37°C. Pick transformants.

Restriction enzymes that do not cut in pJ201:

AarI, AatII, Acc65I, Accl, AclI, AfeI, AgeI, AleI, Alol, ApaBI, AscI, AvrII, BaeI, BamHI, Bbr7I, BbsI, BbvCI, BclI, BclII, BglII, BlnI, BmgBI, BmtI, BplI, BpmI, BsaAI, BsaBI, BsaI, BsaXI, BseRI, BsgI, BsiWI, BsmBI, BsmI, BspEI, BspMI, BsrDI, BsrGI, BstAPI, BstBI, BstEII, BstXI, BstZ17I, Bsu36I, BtgI, ClaI, DraIII, EagI, EarI, EcoICRI, EcoRI, Fall, FseI, FspAI, FspI, HindIII, HpaI, KpnI, MfeI, MluI, MscI, NaeI, NcoI, NdeI, NgoMIV, NheI, NotI, PacI, PfoI, PmeI, PmlI, PpiI, PpuMI, PshAI, PstI, PvuII, RsrII, SacI, SacII, Sall, SanDI, SapI, SbfI, ScaI, Scil, SexAI, SfiI, SgrAI, SnaBI, SpeI, SphI, Sse232I, Sse8647I, StuI, StyI, Swal, TatI, Tth111I, XcmI, XhoI, XmnI, ZraI

Restriction enzymes that do not cut in pJ204:

AarI, AatII, Acc65I, Accl, AgeI, AleI, Alol, ApaBI, AscI, AsiSI, AvrII, BaeI, BamHI, Bbr7I, BbsI, BbvCI, BclI, BglII, BlnI, BmgBI, BmtI, BplI, BpmI, Bpu10I, BsaAI, BsaBI, BsaI, BsaXI, BseRI, BsgI, BsiWI, BsmBI, BsmI, BspMI, BsrDI, BsrGI, BstAPI, BstBI, BstEII, BstXI, BstZ17I, Bsu36I, BtgI, ClaI, DraIII, EagI, EarI, EcoICRI, EcoNI, EcoRI, Fall, FseI, FspAI, HindIII, HpaI, KpnI, MfeI, MluI, MscI, NaeI, NcoI, NdeI, NgoMIV, NheI, NotI, PacI, PflMI, PfoI, PmeI, PmlI, PpuMI, PshAI, PstI, PvuII, RsrII, SacI, SacII, Sall, SanDI, SapI, SbfI, ScaI, SexAI, SfiI, SgrAI, SnaBI, SpeI, SphI, Sse232I, Sse8647I, StuI, StyI, Swal, Tth111I, XcmI, XhoI, ZraI

Licenses**ATUM Catalog Product Purchases**

By purchasing ATUM Catalog Products, the purchaser agrees these products are for single laboratory/company use only. The Catalog Products shall not be distributed to any other party without the written permission of ATUM.